

# Mutector™ Reagents



# CYP2D6

## Genotyping



***Multiplex***

**Multiplexing**

**Sequencing accuracy**

**Results in less than 4 hours**

## Target SNPs

SNP	Variation	SNP	Variation
CYP2D6 * 2	2850C > A	CYP2D6 * 6	1707delT
CYP2D6 * 3	2549delA	CYP2D6 * 9	2615-17delAAG
CYP2D6 * 4	1846G > A	CYP2D6 * 10	100C > T
CYP2D6 * 5	Gene Deletion	CYP2D6 * 17	1023C > T
CYP2D6	4180G>C	CYP2D6 * 41	2988G > A



## FEATURES

### Clear-Cut Results

- Individual SNP is defined by peak's color & size
- Easy to analysis the allele

### Multiplexing Capability

- All SNPs in different exons can be tested in two tubes

### Flexibility

- Easy to add new SNPs without changing protocol

### Simple Protocol for Same Day Results

- Have results in 3-4 hours



## ORDER INFORMATION

Cat. No.	Kit Description	Size
GP11	Mutector™ Reagents for CYP2D6 Genotyping	32 rxn

*For research use only, not for use in diagnostic procedures*



**888.825.6005**

[www.trimgen.com](http://www.trimgen.com)

34 LOVETON CIRCLE, SUITE 210, SPARKS, MD 21152 USA,  
PHONE: 410.472.1100 FAX: 410-472-1303